



Product overview Dental

Cupola Dental

Indication: small gaps
(1-2 missing teeth)
or curved edentulous
areas (maxillary and
mandibular frontal area)



		Before Swelling*			After Swelling**			Swelling time**
		Volume	Projection	Diameter	Volume	Projection	Diameter	
400-2035	Cupola Dental 0.35 ml	0.05 ml	3 mm	6 mm	0.35 ml	5,6 mm	9 mm	40 Days


* without silicone shell

** in vitro in 0.9% NaCl-Sol.

Cylinder Dental

Indication: straight edentulous areas
(lateral mandible and maxilla)



		Before Swelling*			After Swelling**			Swelling time**
		Volume	Length	Diameter	Volume	Length	Diameter	
400-1024	Cylinder Dental 0.24 ml	0.045 ml	7.5 mm	3 mm	0.24 ml	12 mm	6 mm	20 Days
400-1070	Cylinder Dental 0.7 ml	0.15 ml	12 mm	4 mm	0.7 ml	20 mm	7 mm	40 Days
400-1130	Cylinder Dental 1.3 ml	0.25 ml	13 mm	5 mm	1.3 ml	22 mm	9 mm	50 Days
400-1210	Cylinder Dental 2.1 ml	0.42 ml	15 mm	6 mm	2.1 ml	24 mm	10.5 mm	90 Days

* without silicone shell

** in vitro in 0.9% NaCl-Sol.